





The clean and wet cutlery is given directly into the TD 3000L.

The cutlery is spirally transported upwards through a high-grade sterile granulate and is dried and polished. The installed ventilator prevents that granulate comes out of the machine.

The regular use of TD 3000L keeps your cutlery shining.

### The simple operation of TD 3000L

Whether chromium, stainless steel or even high-grade silver – the clean and wet cutlery is given directly from the dishwasher into the cutlery polishing machine TD 3000L.

Via the internal vibrating conveyor it is spirally transported upwards through a high-grade granular substance. This purely natural product dries and polishes your cutlery and keeps it shining– with just a single run. The granulate has to be exchanged every 4 – 8 weeks depending on the frequency the machine is used. The machine operates very quietly due to a special spring system.

Your cutlery will not only pass the inspection by your guests, it even meets the strict EU-Hygienic Guidelines: a special UVC-lamp sterilizes both, the cutlery and the granulate.

The constant use of the TD 3000L reduces the oxidation of cutlery in a remarkable way and makes it look noticeably better. As manufacturer we offer you a 10-years spare parts-guarantee and service.

### Special technical features

- ▶ low priced model for first-time users
- ▶ space saving
- ▶ ventilator for perfect shine



Option – mobile underframe TD 3000L

### EXTRACT FROM OUR PRICE LIST

ART. NUMBER    ART. DESCRIPTION

PU    NET PRICE IN EURO

0017015	Cutlery polishing machine TD 3000L	1	8.," –00,00 AUD
	Capacity:	3.000 cutlery pieces / hour	
	Dimensions:	(l x w x h) approx. 60 x 57 x 40 cm	
0014768	Mobile underframe TD 3000L with extendable drawer	1	
0016462	Mobile underframe TD 3000L without drawer	1	
0016764	Cutlery catcher TD 3000L with sound absorbing effect	1	
0017212	Basket holder TD 3000L	1	
0017376	Original Dörr polishing granulate SG 3000 L	1	
0015204	Cleaning alcohol	1	

03/2011